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
	Phosphor CCD	CMOS	Hybrid
Sensitivity	◎	○	◎
Dynamic Range	△	△	◎
S/N	○	△	◎
Energy resolution	×	×	◎
Data quality	○	○	◎
Active area size	○	◎	△
Multiple module	○	◎	△
Pulse X-ray	integrate	integrate	×
Cost	△	○	○

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	RadEye200	Shadow-Box	RDI	Xineos	Dexela
Chip/Manufacture	Teledyne	Teledyne	Teledyne	Teledyne	Perkin-Elmer
Type		1280HS	RDI	1313-EO/FL	Dexela1512
Area Size(mm)	100x100	128x128	94x148	130x130	115x145
Pixel Size(um)	100x100	100x100	100x100	100x100	75x75
note			Water cooled	FL is 5 times more sensitive than EO. EO has high dynamic range.	No cooling



Tube Voltage: 90kV
Current: 160uA

Cu0.06mm+Al0.5mm

Geometry:
SID:265mm

Rayence
Dexela

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- Rayence 1515CF

Scintillator	Gd2O2S:Tb
Active Area	72.0 x 72.0mm
Active Pixels	3,000 x 2,999px 1,500 x 1,499px
Pixel Pitch	24 / 48 μm
Limiting Resolution	20.8 lp/mm
A/D Conversion	14 bits
Frame Rate	5 / 20 fps
Interface	Camera Link
Energy Range	40~130 Kv
Dimension(W x L x H)	146 x 215 x 29.8 mm
Weight	2.04 kg (4.5lbs)
Certification	CE

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Copyright notice

Copyright © 2013 — Rigaku Corporation and its Global Subsidiaries. All Rights Reserved.

The textual, photographic, video, audio, and combined audiovisual materials and documents resulting from the promotional or educational activities of Rigaku Corporation and its Global Subsidiaries, including the material contained on the various Corporate Web sites, are protected under U.S., Japanese and international laws as copyrighted works. Anyone who displays, reproduces, copies, creates derivative works, or sells our textual, photographic, multimedia, PowerPoint, video or audiovisual programs for commercial or non-commercial purposes without our permission violates the copyright laws and is liable for copyright infringement.

Likewise, Rigaku Corporation (and its Global Subsidiaries) trademarks and service marks are protected by state, federal, and international trademark laws. Any person who uses our marks for commercial or non-commercial purposes without our permission on goods or services in such a way that it dilutes the distinctive quality of our marks or that creates the likelihood of confusion with our marks is liable for trademark infringement.

Liability for copyright or trademark infringement involves the potential for significant civil damages, including, in particular cases, statutory damages, liability for up to three times actual damages, and attorney's fees.

If you have any questions regarding this notice, please send an e-mail to info@rigaku.com.

Copyright confidentiality notice

This presentation is the property of Rigaku Corporation and its subsidiaries (the “Company”) and is strictly confidential. It contains information intended only for the person(s) to whom it is transmitted. With receipt of this information, recipient(s) acknowledges and agrees that: (1) this document is not intended to be distributed, and if distributed inadvertently, will be returned to the Company as soon as possible; (2) the recipient will not copy, fax, reproduce, divulge, or distribute this confidential information, in whole or in part, without the express written consent of the Company; (3) all of the information herein will be treated as confidential material with no less care than that afforded to its own confidential material.